

Test Plan v0.1							
Test ID	Test Case Name						
1	Verify small signal amplification through amplifier stages.						
2	Verify frequency response through filter/amplifier stages.						
3	Verify bias removal from reference probe.						
4	Verify visible muscular and brain waves through scope.						
5	Verify ADC sampling.						
6	Test input impedance of probes using a known working system.						
7	Verify the sound notification component.						
8	Verify the detection of contact by the device.						
9	Verify detection of brain waves, sleep stages, and accumulation of sleep from recorded data.						
10	Verify device records detected brainwaves during a session.						
11	Verify charging of battery via micro usb connection.						
12	Verify case water resistance.						
1b	Verify detection of contact and initial brainwaves after the device is attached to the head.						
2b	Verify detection of sleep stage accumulation and user wakeup.						
3b	Verify device powers down on loss of head contact.						
4b	Verify battrey life requirements.						

<b>Test Author:</b> Julia Filipchuk		<b>Test ID:</b> 1	
<b>Test Case Name:</b> Verify small signal amplification through amplifier stages.		<b>Type:</b> <input checked="" type="checkbox"/> White Box	
<b>Description:</b> Generate a small intensity signal to test the amplification range of the device circuits. We will use a function generator in combination with a voltage divider to generate peak to peak signals of ~1uV and ~20uV. Check amplification is 23 after AD620 stage, 598 (23*26) after stage 1, and 15548 (23*26*26) after stage 2 filter/amplifier. [' 1.0μ', expect, ' 23.0μ', ' 0.598m', ' 15.5m'] [' 20.0μ', expect, ' 0.46m', ' 12.0m', ' 311.0m']		<input type="checkbox"/> Black Box	
<b>Tester Name:</b>		<b>Date:</b>	
<b>HW/SW Version:</b> Breadboard/Solder		<b>Time:</b>	
<b>Test Setup:</b> Device with function generator connected to attenuated input. Function generator outputting a sin wave at 30 Hz. Scope available for signal verification. Power supply for board power.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b> <b>Fail</b> <b>Comments</b>
1	Enable power to board.	Check supplied current is reasonable.	<input type="checkbox"/> <input type="checkbox"/>
2	Connect to attenuated input 1.	N/A	<input type="checkbox"/> <input type="checkbox"/>
3	Measure output of AD620.	Expect ~23μV peak-to-peak signal strength.	<input type="checkbox"/> <input type="checkbox"/>
4	Measure output of stage 1 filter.	Expect ~598μV peak-to-peak signal strength.	<input type="checkbox"/> <input type="checkbox"/>
5	Measure output of stage 2 filter.	Expect ~15.5mV peak-to-peak signal strength.	<input type="checkbox"/> <input type="checkbox"/>
6	Connect to attenuated input 2.	N/A	<input type="checkbox"/> <input type="checkbox"/>
7	Measure output of AD620.	Expect 0.46mV peak-to-peak signal strength.	<input type="checkbox"/> <input type="checkbox"/>
8	Measure output of stage 1 filter.	Expect 12.0mV peak-to-peak signal strength.	<input type="checkbox"/> <input type="checkbox"/>
9	Measure output of stage 2 filter.	Expect 311mV peak-to-peak signal strength.	<input type="checkbox"/> <input type="checkbox"/>
<b>Result:</b>			

  

<b>Test Author:</b> Julia Filipchuk		<b>Test ID:</b> 2	
<b>Test Case Name:</b> Verify frequency response through filter/amplifier stages.		<b>Type:</b> <input checked="" type="checkbox"/> White Box	
<b>Description:</b> Generate a range of input frequencies to test frequency response. High pass cut-off frequency is 0.16Hz. Low-pass cut-off frequency is 47 Hz. Expect 0.5 of Vpp at cut-off frequencies. With two filter stages $1/\sqrt{2} * 1/\sqrt{2} = 0.5$ total attenuation.		<input type="checkbox"/> Black Box	
<b>Tester Name:</b>		<b>Date:</b>	
<b>HW/SW Version:</b> Breadboard/Solder		<b>Time:</b>	
<b>Test Setup:</b> Device with function generator connected to attenuated input 2. Power supply connected to board. Output voltage constant. Scope at output of stage 2 filter.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b> <b>Fail</b> <b>Comments</b>
1	Enable power to board.	Check supplied current is reasonable.	<input type="checkbox"/> <input type="checkbox"/>
2	Connect to attenuated input 2.	N/A	<input type="checkbox"/> <input type="checkbox"/>
3	Generator frequency to 47 Hz.	Output voltage is ~50% of 311mV (156mV).	<input type="checkbox"/> <input type="checkbox"/>
4	Generator frequency to 20 Hz.	Output voltage is ~100% of 311mV.	<input type="checkbox"/> <input type="checkbox"/>
5	Generator frequency to 0.16 Hz.	Output voltage is ~50% of 311mV (156mV).	<input type="checkbox"/> <input type="checkbox"/>
6			
7			
8			
9			
<b>Result:</b>			

	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	3	
	<b>Test Case Name:</b>	Verify bias removal from reference probe.	<b>Type:</b>	<input checked="" type="checkbox"/> White Box	
	<b>Description:</b>	Verify bias is removing noise via the reference probe. Use the function generator a noisy wave and mirror of that same noise. Use the noise as the reference probe.		<input type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>	Device board with external power supply. Probes connected to head with earclip to bias. Output of filter stage 2 to scope. Function generator channel 2 attached to bias.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Disconnect earclip.	View output on scope as baseline.	<input type="checkbox"/>	<input type="checkbox"/>	
2	Connect earclip.	Verify output on scope improves as noise is removed.			
3					
4					
5					
6					
7					
8					
9					
	<b>Result:</b>				
	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	4	
	<b>Test Case Name:</b>	Verify visible muscular and brain waves through scope.	<b>Type:</b>	<input checked="" type="checkbox"/> White Box	
	<b>Description:</b>	Verify visible electromyography (EMG) muscular response and electroencephalogram (EEG) brain waves through amplifier stages.		<input type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>	Device board with external power supply. Probes connected to head with earclip to bias. Output of filter stage 2 to scope.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Verify baseline waves on scope.	Baseline should be int mV range with only low frequencies.	<input type="checkbox"/>	<input type="checkbox"/>	
2	Conduct Forehead scrunch.	Notice EMG signal on the scope.	<input type="checkbox"/>	<input type="checkbox"/>	
3	Conduct Forehead scrunch.	Notice EMG signal on the scope.	<input type="checkbox"/>	<input type="checkbox"/>	
4	Conduct Forehead scrunch.	Notice EMG signal on the scope.	<input type="checkbox"/>	<input type="checkbox"/>	
5					
6					
7					
8					
9					
	<b>Result:</b>				

	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	5	
	<b>Test Case Name:</b>	Verify ADC sampling.	<b>Type:</b>	<input checked="" type="checkbox"/> White Box	
	<b>Description:</b>	We will verify sampling data correctly with the ADC. We will generate an input wave and record the sampled signal. Check the recorded signal has expected values.		<input type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>	Beaglebone board with external power supply. Input into ADC connected to funct			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Generate input sin 0-3V @ 30Hz.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
2	Record 30s.	Recording should have benn created.	<input type="checkbox"/>	<input type="checkbox"/>	
3	Compare recording to input.	Recording should match input.	<input type="checkbox"/>	<input type="checkbox"/>	
4					
5					
6					
7					
8					
9					
	<b>Result:</b>				
	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	6	
	<b>Test Case Name:</b>	Test input impedance of probes using a known working system.	<b>Type:</b>	<input checked="" type="checkbox"/> White Box	
	<b>Description:</b>	Test probes are able to measure data on a working system. We will use the OpenBCI system to verify probes are usable to record brain signals.		<input type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>	Connect probes to OpenBCI system and earprobes. Probes on forehead with earclip attached.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Observe recording on OpeBCI system.	Confirm brainwaves are detectible with probes.	<input type="checkbox"/>	<input type="checkbox"/>	
2					
3					
4					
5					
6					
7					
8					
9					
	<b>Result:</b>				

	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	7	
	<b>Test Case Name:</b>	Verify the sound notification component.	<b>Type:</b>	<input checked="" type="checkbox"/> White Box	
	<b>Description:</b>	Use a test program to generate the range of tones we will use for the device. Confirm all are heard clearly.		<input type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>	Device with speaker circuit connected. External power supply to board.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Start tone test program.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
2	Wait for next tone.	Confirm contact tone plays.	<input type="checkbox"/>	<input type="checkbox"/>	
3	Wait for next tone.	Confirm ready tone plays.	<input type="checkbox"/>	<input type="checkbox"/>	
4	Wait for next tone.	Confirm error tone plays.	<input type="checkbox"/>	<input type="checkbox"/>	
5	Wait for next tone.	Confirm disconnect tone plays.	<input type="checkbox"/>	<input type="checkbox"/>	
6	Wait for next tone.	Confirm wake-up tone plays.	<input type="checkbox"/>	<input type="checkbox"/>	
7					
8					
9					
	<b>Result:</b>				
	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	8	
	<b>Test Case Name:</b>	Verify the detection of contact by the device.	<b>Type:</b>	<input checked="" type="checkbox"/> White Box	
	<b>Description:</b>	Verify detection of contact by capacitive touch to the device. Utilize a test program that plays a tone when contact is detected past a threshold. Test condition with simulated sweat as well as a clean forehead. Test for no false positives in a cloth bag.		<input type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>				
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Turn on device.	No tone is played after 1 minute.	<input type="checkbox"/>	<input type="checkbox"/>	
2	Place device on forehead.	Device makes notification tone within 1 second.	<input type="checkbox"/>	<input type="checkbox"/>	
3	Remove device.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
4	Spray salt water on forehead.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
5	Place device on forehead.	Device makes notification tone within 1 second.	<input type="checkbox"/>	<input type="checkbox"/>	
6	Place device into bag. Shake.	No tone is played after 1 minute.	<input type="checkbox"/>	<input type="checkbox"/>	
7					
8					
9					
	<b>Result:</b>				

	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	9	
	<b>Test Case Name:</b>	Verify detection of brain waves, sleep stages, and accumulation	<b>Type:</b>	<input checked="" type="checkbox"/> White Box	
	<b>Description:</b>	Use a recorded sleep session into our sleep stage detection algorithm. Play into the device and expect sleep stages to match. Verify sleep stages accumulate total sleep and wake is triggered.		<input type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>				
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Run algorithm on recorded file.		<input type="checkbox"/>	<input type="checkbox"/>	
2	Before sleep data.	Brainwaves but not sleep detected. Wake detected.	<input type="checkbox"/>	<input type="checkbox"/>	
3	At sleep.	Stage N1 detected.	<input type="checkbox"/>	<input type="checkbox"/>	
4	At deeper sleep.	Stage N2 detected.	<input type="checkbox"/>	<input type="checkbox"/>	
5	Later.	Accumulation of N1, N2 stages trigger wake-up.	<input type="checkbox"/>	<input type="checkbox"/>	
6					
7					
8					
9					
	<b>Result:</b>				
	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	10	
	<b>Test Case Name:</b>	Verify device records detected brainwaves during a session.	<b>Type:</b>	<input checked="" type="checkbox"/> White Box	
	<b>Description:</b>	Test if data recording of detected brainwaves functions.		<input type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>	Device must have room on SD card.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Connect device to head.	Hear signal device is connected.	<input type="checkbox"/>	<input type="checkbox"/>	
2	Record 40 seconds.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
3	Scrunch forehead 10 times.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
4	Record 40 seconds.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
5	Stop recording.		<input type="checkbox"/>	<input type="checkbox"/>	
6	View recored file.	File should exist on SD card.	<input type="checkbox"/>	<input type="checkbox"/>	
7	View data.	Scrunches should be visible as spikes.	<input type="checkbox"/>	<input type="checkbox"/>	
8					
9					
	<b>Result:</b>				

	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	11	
	<b>Test Case Name:</b>	Verify charging of battery via micro usb connection.	<b>Type:</b>	<input checked="" type="checkbox"/> White Box	
	<b>Description:</b>	Test if battery can be charged in the system.		<input type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>	Device.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Discharge battery.	Measured voltage under rated 3.6 v indicating ~10% charge.	<input type="checkbox"/>	<input type="checkbox"/>	
2	Connect to power.	Battery starts charging.	<input type="checkbox"/>	<input type="checkbox"/>	
3	Wait for 1h.	No fire. Device is still on.	<input type="checkbox"/>	<input type="checkbox"/>	
4	Disconnect power.	Device still on.	<input type="checkbox"/>	<input type="checkbox"/>	
5	Measure battery voltage.	Measured voltage above 3.7 V indicating > 60% charge.	<input type="checkbox"/>	<input type="checkbox"/>	
6					
7					
8					
9					
	<b>Result:</b>				
	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	12	
	<b>Test Case Name:</b>	Verify case water resistance.	<b>Type:</b>	<input checked="" type="checkbox"/> White Box	
	<b>Description:</b>	Place in a moist tupperware with a moisture test strip inside. Spray with a spray bottle 3x. Leave overnight.		<input type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>	Case with out main board but with all external connections.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Place Test Strip Inside.	Test strip is unactivated (dry).	<input type="checkbox"/>	<input type="checkbox"/>	
2	Close case.	No gaps on case.	<input type="checkbox"/>	<input type="checkbox"/>	
3	Spray 3x with water.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
4	Close in tupperware container.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
5	Wait 24 hours.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
6	Check if test strip has been activated.	Test strip should signal that himidity was low.	<input type="checkbox"/>	<input type="checkbox"/>	
7					
8					
9					
	<b>Result:</b>				

	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	1b	
	<b>Test Case Name:</b>	Verify detection of contact and initial brainwaves after the device	<b>Type:</b>	<input type="checkbox"/> White Box	
	<b>Description:</b>	In order to test correct device positioning we will have a user place the device on the forehead. Audible confirmation if in contact and brainwaves in threshold. Audible alert if contact but the brainwaves are undetectable or under threshold.		<input checked="" type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>	Charged device with clean probes.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Place on head.	Listen for contact confirmation tone within 10 seconds.	<input type="checkbox"/>	<input type="checkbox"/>	
2	Wait.	Listen for brainwave confirmation tone within 40 seconds.	<input type="checkbox"/>	<input type="checkbox"/>	
3					
4					
5					
6					
7					
8					
9					
	<b>Result:</b>				
	<b>Test Author:</b>	Julia Filipchuk	<b>Test ID:</b>	2b	
	<b>Test Case Name:</b>	Verify detection of sleep stage accumulation and user wakeup.	<b>Type:</b>	<input type="checkbox"/> White Box	
	<b>Description:</b>	Verify core functionality of device. User will take a nap with the expectation that device will take them.		<input checked="" type="checkbox"/> Black Box	
	<b>Tester Name:</b>		<b>Date:</b>		
	<b>HW/SW Version:</b>	Breadboard/Solder	<b>Time:</b>		
	<b>Test Setup:</b>	Charged device with clean probes. Second person to wake user. Tired main device user.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Place device on forehead.	Contact tone.	<input type="checkbox"/>	<input type="checkbox"/>	
2	Device is ready.	Confirmation tone.	<input type="checkbox"/>	<input type="checkbox"/>	
3	User goes to sleep.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
4	Wait ~20+ minutes.	Device plays wakeup tone after ~10-20 minutes asleep.	<input type="checkbox"/>	<input type="checkbox"/>	
5					
6					
7					
8					
9					
	<b>Result:</b>				



		<b>Test Author:</b> Julia Filipchuk	<b>Test ID:</b> 3b		
		<b>Test Case Name:</b> Verify device powers down on loss of head contact.	<b>Type:</b> <input type="checkbox"/> White Box		
		<b>Description:</b> Check the device powers off after removing from head contact.	<input checked="" type="checkbox"/> Black Box		
		<b>Tester Name:</b>	<b>Date:</b>		
		<b>HW/SW Version:</b> Breadboard/Solder	<b>Time:</b>		
		<b>Test Setup:</b> Connect device to head and allow it to start recording.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Remove device from head.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
2	Verify power is off.	Power off within 30 seconds.	<input type="checkbox"/>	<input type="checkbox"/>	
3					
4					
5					
6					
7					
8					
9					
		<b>Result:</b>			
		<b>Test Author:</b> Julia Filipchuk	<b>Test ID:</b> 4b		
		<b>Test Case Name:</b> Verify battrey life requirements.	<b>Type:</b> <input type="checkbox"/> White Box		
		<b>Description:</b> Run device on battery for 3 days. Using it for wake-up 6 times.	<input checked="" type="checkbox"/> Black Box		
		<b>Tester Name:</b>	<b>Date:</b>		
		<b>HW/SW Version:</b> Manufactured Board in Case	<b>Time:</b>		
		<b>Test Setup:</b> Device in case with strap. Charge device fully. Clean electrodes.			
<b>Step</b>	<b>Action</b>	<b>Expected Result</b>	<b>Pass</b>	<b>Fail</b>	<b>Comments</b>
1	Wait until sleep.	N/A	<input type="checkbox"/>	<input type="checkbox"/>	
2	Place on forehead.	Device beeps within 30 seconds indicating ready.	<input type="checkbox"/>	<input type="checkbox"/>	
3	Fall asleep.	Device wakes you up in ~20 minutes from you falling asleep.	<input type="checkbox"/>	<input type="checkbox"/>	
4	Place device in bag.	Device is silent.	<input type="checkbox"/>	<input type="checkbox"/>	
5	Wait ~12 hours. Repeat 1-5.		<input type="checkbox"/>	<input type="checkbox"/>	
6	Wait ~12 hours. Repeat 1-5.		<input type="checkbox"/>	<input type="checkbox"/>	
7	Wait ~12 hours. Repeat 1-5.		<input type="checkbox"/>	<input type="checkbox"/>	
8	Wait ~12 hours. Repeat 1-5.		<input type="checkbox"/>	<input type="checkbox"/>	
9	Wait ~12 hours. Repeat 1-5.	Device still functions after 6 uses in 3 days.	<input type="checkbox"/>	<input type="checkbox"/>	
		<b>Result:</b>			